EN

Carefully read the following instructions before using the guidewire.

Description
The RocaWire guidewires are sterile and are intended for single-use only. Please refer to the packaging for the specifications of each guidewire.

Indications
RocaWire guidewires are medical devices designed to facilitate the positioning of endourological devices during endourological procedures. The Nephro RocaWire guidewires can also be used to facilitate the positioning of devices during percutaneous nephrostomy procedures.

Contraindications
None known of at the present time.

Possible complications
Complications associated with surgical intervention. Complications associated with any endourological procedure. E.g.: haemorrhaging, lesions/perforation of tissue, perforation of urinary tracts, infection, dislodgement of catheter or stent, vesicoureteral reflux, urinary retention.

- Precautions:**
The following usage precautions should be carefully observed to ensure the proper function of the device and optimal comfort and safety for the patient:
- Do not use the guidewire if the packaging is damaged as this could compromise guidewire integrity and/or sterility.
 - Do not use the guidewire if it has reached or passed its expiry date.
 - Verify the integrity of the guidewire before use. Do not use the guidewire if it is damaged either before or during use.
 - The user should be thoroughly trained in endoscopic techniques before using this device. Any incorrect use, kinking, excessive bending, excessive progression speed, deterioration of the external coating, the distal tip, or any forms of use other than those indicated in the usage instructions or that may be contrary to these precautions, may lead to the deterioration of the guidewire.
 - Hydrophilic guidewires must not be inserted into a metal tube, a needle or any metal object.
 - The use of the PTFE-coated guidewires with metal instruments is only permissible if precautions are taken to prevent damage.
 - During simultaneous use of the guidewire and other products, medical devices and instruments, the surgeon must verify the characteristics and compatibility of all devices to avoid damage to the guidewire or to other devices.
 - When used simultaneously with other devices that could damage the guidewire, the practitioner must ensure the condition of the guidewire throughout the entire radiocopy or endoscopy procedure.
 - The guidewires must be manipulated and used with care, by means of radiocopy or endoscopic monitoring, so as to prevent damage to the guidewire and to ensure the perfectly atraumatic progression of the guidewire towards the target area. If there is any doubt whatsoever, the practitioner must withdraw the guidewire under radiocopy and/or endoscopy, with the extreme care.
 - Pay special attention during use of the device in combination with a laser fibre: the laser could easily damage the device.
 - If the guidewire appears to be damaged during the procedure, withdraw it immediately whilst taking the necessary care to prevent damage. Ensure that no part of the guidewire remains inside the body of the patient.
 - Antiseptic solutions, in particular those containing alcohol, must be avoided as they may alter the surface of the guidewire.
 - Use of guidewire forceps is not recommended. It should be noted that RocaWire guidewires have a metal core; use with inappropriate instruments is not recommended.
 - The guidewires must be used under antiseptic conditions, in the operating theatre, in compliance with the regulations in force in the respective county.

Instructions for use

Hydrophilic models (references ROWP / ROWR):

Withdraw the RocaWire Platinum guidewire from its peel-open pouch. Using a syringe, inject a physiological solution into the protective sheath to reveal the hydrophilic coating. Then, when moistened, remove the RocaWire guidewire from its protective sheath without forcing. If it does not come out easily, inject the physiological solution a second time and wait a few seconds before withdrawing the guidewire from its protective sheath.

With prolonged usage, the surface of the guidewire may become dry and less slippery. If this is the case, moisten the guidewire again using the physiological solution. Due to the very slippery nature of the hydrophilic models, the guidewire must be held firmly in place during use, by means of a torque device if necessary.

Inspect the guidewire before use. Insert the guidewire into the patient slowly, whilst monitoring it visually or by means of radiocopy in order to prevent any lesions or damage.

Other models (references ROWG / ROWS):

Withdraw the RocaWire guidewire from its peel-open pouch. Injection of the physiological solution is not necessary. Inspect the guidewire before use. Insert the guidewire into the patient slowly, whilst monitoring it visually or by means of radiocopy in order to prevent any lesions or damage.

Usage and storage conditions:

The set should be stored and transported in a normal environment, i.e. away from extreme temperatures and humidity. Do not use if the sterile packaging is damaged and/or deteriorated. Single-use device: do not reuse, reprocess or resterilise. This could compromise the integrity of the device and/or cause malfunctions that could potentially lead to injury, illness or death of the patient. In addition, there is a risk of contaminating the device, which could cause infections in the patient or cross-infection with another patient, thus leading to injury, illness or death. The manufacturer will not be held liable in such events. After use, discard the device and packaging in compliance with hospital, local authority and/or government policies.

Liabilities and Guarantee:

PROMEPLA and its affiliates declare to the first purchaser of this device that reasonable care has been taken in the design and manufacture thereof. If damage is suspected, please contact PROMEPLA and its representatives. The present limited guarantee dispenses all other guarantees, express or implied, including but not limited to any implied guarantees concerning merchantability or suitability for specific purposes. The liability of PROMEPLA and its companies affiliated by virtue of this limited product guarantee shall be limited to the replacement of defective products. Under no circumstances shall **PROMEPLA and its affiliates be held liable for direct, indirect, incidental and/or consequential damages of any kind arising from the use and/or handling of this device.**

 	 	 	
		Manufacturer : ROCAMED France – ZI de Signes Allée de Stockholm - 83870 Signes – FRANCE Tel. +33 494 902 100 - Fax. +33 494 986 055	
		This device is not made with natural rubber latex	
 	 	 	CE marked on: 20 September 2013

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