## EN

### Description

The RocaTub catheters are sterile and for single use only. Please refer to labeling for the specifications of each product.

ROCAMED's RocaTub ureteral catheter series offers a wide range of catheters with various distal tip shapes.

### Indications

RocaTub catheters are intended for use during ureteral catheterisation, drainage, opacification of the upper urinary tract (retrograde ureteropyelography) and flushing procedures.

#### Contraindications

Any contraindications related to ureteral catheterisation known of at the present time.

#### Possible complications

A number of adverse events have been reported in the literature with indwelling catheters: infection, encrustation, obstruction, migration, vesical irritation, pain and erosion. Events such as ureteral perforation can occur during intraoperative procedures, the most often related to concomitant use of a guidewire.

#### Precautions

The operator should observe the following precautions conscientiously to ensure patient safety and proper functioning of the device.

- Do not use the catheter if the packaging is damaged because this could compromise device integrity and/or sterility.
- Do not use the catheter after its used by date.
- Do not use the catheter if damaged or kinked either before or during use.
- Care should be taken to ensure that the hydrophilic coating is properly activated upon insertion, especially on account of ventilation systems and air conditioning in operating rooms which can rapidly dry the coating. For indwelling ureteral catheters, regular follow-up is necessary after implantation to ensure that the device is working properly.
- For indwelling ureteral catheters, the patient should be informed about the necessary regular follow-up. The patient should be informed of the importance of contacting the physician immediately in the event of dysfunction of the device.
- If infection, encrustation, decrease and/or interruption of drainage is suspected, the catheter must be removed.

#### Instructions for use

Open both the outside and inner sterile pouches and take the RocaTub out using usual aseptic techniques. Visually inspect the device for defects and check that the stylet is mobile to ensure that subsequent removal is possible. Likewise, make sure that the stylet does not protrude from the distal tip before catheter insertion. Immerse the distal tip of the catheter in physiological saline or sterile water for a few seconds to activate the hydrophilic coating. **Note:** ventilation systems and air conditioning in operating rooms can rapidly dry the hydrophilic coating. Make sure the hydrophilic coating is properly activated before insertion and re-immers the catheter if necessary. Insert the catheter as per the usual procedure. Catheter insertion under endoscopy is usually monitored visually and/or under fluoroscopy and/or radioscopy. Once introduced and correctly positioned, remove the stylet. The stylet may be removed before insertion if the catheter is passed over a guidewire.

**Warning:** For intraoperative catheters and/or indwelling catheters, catheter features (size, tip, number of apertures, etc.) should be selected to match the shape and size of the ureter, the size of the patient and the desired use. In all cases, it is the surgeon's responsibility to make these choices.

**Injection of contrast medium:** Dilute the contrast media in physiological saline. Immediately inject using the syringe and the injection tip provided with the device. Warning: The injection of contrast medium can result in adverse effects. Care should be taken to follow all the safety recommendations described on the package insert.

**Removal of the device** The device is removed by simply pulling on the proximal end.

### Conditions of use and storage

The device should be stored and transported in a normal environment, i.e. away from extreme temperatures and humidity. Do not use if the sterile packaging is damaged and/or open. For single use only; resterilization of this device jeopardizes sterility, can lead to product dysfunction and cross-contamination. The manufacturer will not be held liable in the case of re-use. Destroy after use according to the regulations in force regarding disposal of potentially infectious waste material.

#### Liability and warranty

Rocamed and affiliates declare to the first purchaser, that reasonable care has been taken when designing and manufacturing this product. If damage is suspected, please contact Rocamed and affiliates. The present limited warranty dispenses all other warranties, express or implied, including but not limited to any implied warranties concerning usability or suitability for specific purposes. This limited guarantee ensures that the responsibility of Rocamed and affiliated companies will be limited to the replacement of defective products. In no instance will Rocamed and affiliates be held liable for direct, indirect, incidental and/or consequential damages of any kind arising from the use and/or handling of the device.

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